




FRED easy



The Automated External Defibrillator  
that accompanies all first aid rescuers in  
their life-saving mission



**SCHILLER**

The Art of Diagnostics

# FRED easy

The FRED easy Automated External Defibrillator is the best life-saving tool for all emergency services: paramedics, ambulances, firefighters....

It is light, rugged and user friendly and therefore adapted to all situations.

FRED easy features the unique pulsed biphasic waveform MULTIPULSE BIOWAVE with patient impedance compensation. This waveform is efficient for adults and children thanks to the automatic recognition of children pads and energy adaptation.

Adapted voice and displayed instructions for first aid rescuers together with a fast rhythm analysis and shock delivery guide the rescuer during the resuscitation process and minimise the hands-off time.

In order to follow the patient throughout the resuscitation process, the intervention data is recorded and can be reviewed on a PC and a report can be generated.

FRED easy is available in a range of versions and options to fulfil the rescuer's needs:

- Semi or fully automatic
- Metronome: providing the correct rhythm for chest compressions
- ECG display
- Manual override: doctors can switch off the AED mode and decide whether they want to defibrillate or not




FRED easy with manual override



FRED easy fully automatic



ECG display



## FRED easy Online

The first defibrillator featuring automatic data transmission. Intervention data and the device configuration can be sent to a server via an Ethernet network.

### Advantages for everyone:

Transmission is easy, fast and safe. The data is stored centrally for direct and immediate access by authorised persons from any computer connected to the network. Multi-user access to the data via any computer connected to the network (hospital, treating physician, cardiologist, rescue services, administration). This system considerably simplifies the handling of medical and administrative patient data (no paper, no memory card).

### LifeDataNet:

The transmitted data is stored on a central server, LifeDataBox, with the LifeDataNet application.

### Two types of data can be transmitted:

- Patient data (events, ECG curve, patient impedance)
- Configuration and maintenance information

### The LifeDataNet application is a simple and user-friendly interface enabling you to:

- see which device has transmitted the data (serial number, identification, transmission date and time)
- download and analyse patient data
- view the device settings (mode, energy settings, options)
- view the maintenance status (self-test result, battery status, software version).

When a first aid vehicle comes back to the station, simply connect the FRED easy to the network and the data is automatically downloaded and the device ready for a new intervention.

Schiller also offers ARGUS PRO LifeCare 2 and DG4000 to perform advanced patient resuscitation and monitoring.



ARGUS PRO LifeCare 2



DEFIGARD 4000



**Asia**  
**SCHILLER Asia-Pacific / Malaysia**  
52200 Kuala Lumpur, Malaysia  
Phone +603 6272 3033  
Fax +603 6272 2030  
sales@schiller.com.my  
www.schiller-asia.com



**Austria**  
**SCHILLER Handelsgesellschaft m.b.H.**  
A-4040 Linz  
Phone +43 732 709 90  
Fax +43 732 757 000  
sales@schiller.at  
www.schiller.at



**China**  
**SCHILLER Int'l Trading (Shanghai) Co. Ltd.**  
200335 Shanghai, China  
Phone +86-21-62099627  
Fax + 86-21-62099623  
sales@schiller.com.my  
www.schiller.cn



**France**  
**SCHILLER Médical S.A.S.**  
F-67162 Wissembourg/Cedex  
Phone +33 3 88 63 36 00  
Fax +33 3 88 94 12 82  
info@schiller.fr  
www.schiller-medical.com



**France (distribution France)**  
**SCHILLER France S.A.S.**  
F-77600 Bussy St Georges  
Phone +33 1 64 66 50 00  
Fax +33 1 64 66 50 10  
infoschiller@schiller-france.fr  
www.schiller-france.com



**Germany**  
**SCHILLER Medizintechnik GmbH**  
D-85521 Ottobrunn  
Phone +49 89 62 99 81-0  
Fax +49 89 609 50 90  
info@schillermed.de  
www.schillermed.de



**Hungary**  
**SCHILLER Diamed Ltd.**  
H-1141 Budapest  
Phone +36 (1) 383-4780 / 460-9491  
Fax +36 (1) 383-4778  
sales@schiller.at  
www.schiller-hungary.hu



**India**  
**SCHILLER Healthcare India Pvt. Ltd.**  
Mumbai - 400 001, India  
Phone +91 22 6152 3333/ 2920 9141  
Fax +91 22 2920 9142  
sales@schillerindia.com  
www.schillerindia.com



**Japan**  
**SCHILLER Japan, Ltd.**  
Hiroshima 734-8551  
Phone +81 82 250 2055  
Fax +81 82 253 1713  
koji.maekawa@schiller.jp  
www.schiller.jp



**Poland**  
**SCHILLER Poland Sp. z o.o.**  
PL-02-729 Warszawa  
Phone +48 22 8432089  
Fax +48 22 8432089  
schiller@schiller.pl  
www.schiller.pl



**Russia & C.I.S.**  
**SCHILLER AG Rep. office**  
125124 Moscow, Russia  
Phone +7 (495) 970 11 33  
Fax +7 (495) 970 11 33  
mail@schiller-ag.com  
www.schiller-cis.com



**Spain**  
**SCHILLER ESPAÑA, S.A.**  
E-28230-Las Rozas/Madrid  
Phone +34 91 713 01 76  
Fax +34 91 355 79 33  
schiller@schiller.es  
www.schiller.es



**Switzerland**  
**SCHILLER-Reomed AG**  
CH-8953 Diefikon  
Phone +41 44 744 30 00  
Fax +41 44 740 37 10  
sales@schiller-reomed.ch  
www.schiller-reomed.ch



**Turkey**  
**SCHILLER TÜRKİYE**  
Okmeydanı-Sisli - İstanbul  
Phone +90 212 210 8681 (pbx)  
Fax +90 212 210 8684  
sales@schiller-turkiye.com  
www.schiller-turkiye.com



**USA**  
**SCHILLER America Inc.**  
Doral, Florida 33172  
Phone +1 786 845 0620  
Fax +1 786 845 06 02  
sales@schilleramerica.com  
www.schilleramerica.com

11/11



Part no.: 2.500378



# SCHILLER

The Art of Diagnostics