



CARDIOVIT MS-2015


Defining the new horizon
of electrocardiography

12-lead high-end ECG system
for resting ECGs

Display
15" full colour touch screen

Speed
Instant-on and ECG printout in less than
3 seconds

DICOM
High-end electrocardiograph
to support DICOM standard



SCHILLER

The Art of Diagnostics

CARDIOVIT MS-2015

The right touch

SCHILLER's MS-2015 is a state-of-the-art electrocardiograph that simplifies diagnostic ECG interpretation and streamlines workflows by delivering high-quality clinical reports when and where you need them.

Prepare to be wowed at your first encounter with the MS-2015 and impressed every time you use it. The MS-2015 is designed for constant use. No matter how hectic the clinical environment – in large hospitals where ECG volume is high and speed and quality are paramount – the MS-2015 is ready to help you to simplify cardiac patient care.

Simply touch the large (15"), high-resolution colour display and you're in business – recording, selecting and printing the highest quality ECGs faster than you can imagine. MS-2015 exports ECG recordings in XML, PDF*¹ and DICOM*¹ format, so you can transmit to any system that supports these industry standards. From accurate and complete data entry to comprehensive data fields that prompt users for required information, every detail of the MS-2015 is designed to streamline ECG testing.

And when you're responsible for hundreds of ECGs every day, attention to detail makes all the difference.

SCHILLER AG
pioneering company
since 1974



Advanced processing... no waiting

- Instant-on start-up provides immediate availability
- ECG processed, interpreted and printed in less than 3 seconds after acquisition

Synchronize time

Auto-set the MS-2015 time with your hospital time master to obtain accurate documentation of your patient's clinical history.

Fast entry

Enter patient data quickly and accurately using the touch screen or by scanning a bar code. You can also download*² patient data and worklist requests from SEMA or HIS.

Clinical excellence

- **SCHILLER Interpretation software "C" for adult and paediatric patients** – one of the best documented algorithms for simultaneous ECG analysis of 12 leads, offering quality and reliability without compromise. Your most valuable tool to get a clinically reliable second opinion.
- **Thrombolysis software** – this option creates a numerical assessment

of the patient's age, gender, main symptoms as well as the ECG measurement, based on which you can determine the probability of acute cardiac ischaemia. This optional program gives you important additional diagnostic information for patients with chest pain.

Never miss a beat*²

Capture and store 4 minutes of uninterrupted 12-lead data to document intermittent arrhythmias.

High-volume storage

- Storage of up to 350 ECGs
- Full-disclosure storage of up to 4 minutes*²



THE ART OF CONNECTIVITY

ECGs that meet your high standards

Not only is the MS-2015 capable of seamless connection to the SEMA-200 information systems, it also adapts to your specific needs and can be expanded as those needs grow. Full connectivity enables you to streamline your ECG workflow and improve efficiency.

- **EMR ready** – The MS-2015 connects easily to your EMR system.
- **Seamless** – Connectivity with SEMA-200 can help clinicians to increase departmental productivity and reduce errors.
- **Simplification** – Automatic storage, transmission and deletion of data.
- **Barcode scanning**^{*2} – The barcode reader option helps reduce errors by automating the input of patient data.

Maintain security

The MS-2015 offers you safe wireless connections via standard WLAN protocols such as 802.11 g, WPA, WPA II or WEP. WLAN communication is also available with invisible SSID.

Spirometry

ECG and spirometry in one unit – this is unique world-wide! The CARDIOVIT MS-2015 can be transformed into a proven table-top spirometer for the measurement, recording and assessment of flow-volume and volume-time curves and the corresponding parameters. Several expiratory and inspiratory tests can be performed and compared with country-specific normal values. Pre/post measurements and interpretation are standard.

^{*1} Export of 12 lead Resting ECG recordings in standard DICOM Waveform Object format to PACS systems and Pdf export is available as part of the SCHILLER Communication Server (SCS).

^{*2} development in progress



Asia
SCHILLER Asia-Pacific / Malaysia
52200 Kuala Lumpur, Malaysia
Phone +603 6272 3033
Fax +603 6272 2030
sales@schiller.com.my
www.schiller-asia.com



Austria
SCHILLER Handelsgesellschaft m.b.H.
A-4040 Linz
Phone +43 732 709 90
Fax +43 732 757 000
sales@schiller.at
www.schiller.at



China
SCHILLER Int'l Trading (Shanghai) Co. Ltd.
200335 Shanghai, China
Phone +86-21-62099627
Fax + 86-21-62099623
sales@schiller.com.my
www.schiller.cn



France
SCHILLER Médical S.A.S.
F-67162 Wissembourg/Cedex
Phone +33 3 88 63 36 00
Fax +33 3 88 94 12 82
info@schiller.fr
www.schiller-medical.com



France (distribution France)
SCHILLER France S.A.S.
F-77600 Bussy St Georges
Phone +33 1 64 66 50 00
Fax +33 1 64 66 50 10
infoschiller@schiller-france.fr
www.schiller-france.com



Germany
SCHILLER Medizintechnik GmbH
D-85521 Ottobrunn
Phone +49 89 62 99 81-0
Fax +49 89 609 50 90
info@schillermed.de
www.schillermed.de



Hungary
SCHILLER Diamed Ltd.
H-1141 Budapest
Phone +36 (1) 383-4780 / 460-9491
Fax +36 (1) 383-4778
sales@schiller.at
www.schiller-hungary.hu



India
SCHILLER Healthcare India Pvt. Ltd.
Mumbai - 400 001, India
Phone +91 22 6152 3333 / 2826 3520
Fax +91 22 2826 3525
sales@schillerindia.com
www.schillerindia.com



Japan
SCHILLER Japan, Ltd.
Hiroshima 734-8551
Phone +81 82 250 2055
Fax +81 82 253 1713
koji.maekawa@schiller.jp
www.schiller.jp



Poland
SCHILLER Poland Sp. z o.o.
PL-02-729 Warszawa
Phone +48 22 8432089
Fax +48 22 8432089
schiller@schiller.pl
www.schiller.pl



Russia & C.I.S.
SCHILLER AG Rep. office
125124 Moscow, Russia
Phone +7 (495) 970 11 33
Fax +7 (495) 970 11 33
mail@schiller-ag.com
www.schiller-cis.com



Spain
SCHILLER ESPAÑA, S.A.
E-28230-Las Rozas/Madrid
Phone +34 91 713 01 76
Fax +34 91 355 79 33
schiller@schiller.es
www.schiller.es




Switzerland
SCHILLER-Reomed AG
CH-8953 Dietikon
Phone +41 44 744 30 00
Fax +41 44 740 37 10
sales@schiller-reomed.ch
www.schiller-reomed.ch



Turkey
SCHILLER TÜRKİYE
Okmeydani-Sisli - Istanbul
Phone +90 212 210 8681 (pbx)
Fax +90 212 210 8684
sales@schiller-turkiye.com
www.schiller-turkiye.com



USA
SCHILLER America Inc.
Doral, Florida 33122
Phone +1 786 845 0620
Fax +1 786 845 06 02
sales@schilleramerica.com
www.schilleramerica.com



SCHILLER

The Art of Diagnostics